

3. BOND (cont.)

opal[®] BOND™ FLOW

FLOWABLE LIGHT-CURED ADHESIVE

- Features optimal flow characteristics that ease delivery and placement and help prevent slumping
- Ideal for permanent retainers, bite ramps, occlusal buildups, and indirect bonding
- Available in two shades: A2 and Blue
- Compatible with any light-cured bonding system
- Kind to opposing enamel and has a smooth finish



opal[®] BAND™

ORTHODONTIC BAND CEMENT

- Triple-cure formula provides the ideal bond strength for multiple uses
- Dual-barrel syringe delivery allows for precise placement and the convenience of no mixing
- Glass ionomer fillers provide a better cure, as well as the added benefit of fluoride release
- Ideal viscosity allows the material to flow easily through the mixing tip and gives complete control when applying
- Bright blue-colored cement allows for easy removal and cleanup



888.863.5883 OPALORTHODONTICS.COM

© 2014 Ultradent Products, Inc. All rights reserved.

95653.4 060914



— BONDING ESSENTIALS FROM OPAL ORTHODONTICS —

COMPLETELY SIMPLE

ETCH • PRIME • BOND

With its complete line of adhesive products, Opal Orthodontics offers a simple, three-step bonding system that ensures consistent, reliable results.



opal[®]
ADHESIVES

1. ETCH

opal[®] ETCH[™]



35% PHOSPHORIC ACID ETCHING SOLUTION

- Brilliant blue color ensures safe, accurate placement
- Optimal viscosity penetrates fissures without migrating
- Syringe and tip offer controlled delivery, protecting soft tissues
- Proprietary surfactant rinses cleanly and completely

2. PRIME

opal[®] SEAL[™]



FLUORIDE RELEASING AND RECHARGING PRIMER

- Releases and recharges fluoride
- 38% filled with substantial glass ionomer plus nanofillers for long-lasting strength
- Nonyellowing and stain resistant
- Drying agent seeks out moisture and draws resin in, ensuring superior bonding and fluoride uptake
- Fluorescent properties ease reapplication and removal

3. BOND

opal[®] BOND[™] MV



MEDIUM-VISCOSITY LIGHT-CURED BONDING ADHESIVE

- Easy to butter into bracket mesh yet tacky enough to minimize drift
- Can be delivered in compules or directly via Opal[®] Orthodontics' unique, economical contra-angle syringe
- Proprietary loading process creates virtually no "run-on"
- Visibility under a black light facilitates easy removal
- Extended working time
- Easy flash cleanup

opal[®] BOND[™]

VISCOUS LIGHT-CURED BONDING ADHESIVE

- High viscosity allows easy positioning and flash removal
- Super tacky, viscous formula helps eliminate bracket drift while providing ample working time
- Convenient, familiar compule delivery

